

## MSCprime<sup>®</sup>, Primary MSC Supplement

### Product Overview:

*Mesenchymal stem cells (MSCs) can be expanded in serum-free media, which is a critical step towards clinical application. However, isolation of the primary MSCs from human tissues in serum-free media suffers from low efficiency due to the lack of cell adhesion proteins.*

**MSCprime<sup>®</sup>** is developed for the isolation of primary MSC from human tissues when used with StemRD's **MesenGro<sup>®</sup>** (now available as **StemGro<sup>®</sup>** hMSC Medium from Corning) a serum-free and chemically-defined MSC medium (#40-410-KIT, Corning-Mediatech). Its features include:

1. Increasing plating efficiency of primary MSCs
2. Resulting in more homogeneous and morphologically-desirable MSCs
3. Serum-free (no human or animal serum added)
4. Chemically-defined (all ingredients are either purified or recombinant)
5. Better lot-to-lot consistency than serum or other undefined supplements

### Package size and storage:

**MSCprime<sup>®</sup>** is a concentrated solution (25 X) and shipped in three package sizes frozen:

	<b>Catalog #</b>	<b>Size (25 X)</b>	<b>Diluted Volume</b>
<b>MSCprime<sup>®</sup></b>	MP-001	1 mL	25 mL
	MP-005	5 mL	125 mL
	MP-025	25 mL	625 mL

- **Storage and Shelf life:** **MSCprime<sup>®</sup>** can be stored at -80°C for 6 months. Diluted **MSCprime<sup>®</sup>** can be stored at 2 to 8°C for 2 weeks.
- **Avoid freeze and thaw.**

### Additional key reagents/supplies needed:

- **MesenGro<sup>®</sup>**, serum-free and chemically-defined MSC medium (**StemRD** cat# MGro-500)
- **CellBIND<sup>®</sup>** tissue culture flasks or dishes (e.g., T25 flasks, **Corning** cat# 3289 )


**FOR RESEARCH USE ONLY. NOT FOR HUMAN THERAPEUTIC OR DIAGNOSTIC USES.**

## Application Instruction:

1. Prior to use, thaw **MSCprime**<sup>®</sup> at 2 to 8°C overnight.
2. Add primary cell isolate containing MSCs (e.g., MNC layer of bone marrow aspirate) into 10X volume of complete **MesenGro**<sup>®</sup> serum-free medium, spin cells down at 1500 rpm for 10 min.
3. During cell spinning, dilute **MSCprime**<sup>®</sup> into complete **MesenGro**<sup>®</sup> medium at a ratio of **1:24**.
4. Resuspend the cell pellet in the diluted **MSCprime**<sup>®</sup>/**MesenGro**<sup>®</sup> medium and transfer into a **CellBIND**<sup>®</sup> tissue culture flask. Place the flask in 37°C, 5% CO<sub>2</sub> incubator.
5. Change medium to complete **MesenGro**<sup>®</sup> in 2 to 6 days. The optimal time of the medium change to **MesenGro**<sup>®</sup> may vary between tissue types. We have successfully isolated MSCs from bone marrow after 2 days and umbilical cord after 6 days in **MSCprime**<sup>®</sup>/**MesenGro**<sup>®</sup>. The appearance of MSC-like colonies can be used as a guide.


## Expected Results:

1. Isolation of MSC-like cells from human umbilical cord by using **MSCprime**<sup>®</sup>:



*Fig 1. Colonies of MSC-like cells isolated with **MSCprime**<sup>®</sup>/**MesenGro**<sup>®</sup> medium are denser, while the cells are more homogeneous and smaller than with FBS.*

2. **MSCprime**<sup>®</sup>-isolated homogeneous/smaller MSCs subsequently expanded in **MesenGro**<sup>®</sup>:



*Fig. 2. MSCs isolated by **MSCprime**<sup>®</sup> expand faster in **MesenGro**<sup>®</sup> with more desirable morphology than in FBS.*

**FOR RESEARCH USE ONLY. NOT FOR HUMAN THERAPEUTIC OR DIAGNOSTIC USES.**